



PEMF Prime Academy

Certified PEMF Expert


Module 4

Understanding and Implementing Heart Rate Variability (HRV)

PEMF Prime Academy

Heart Rate Variability

- ◆ Heart Rate (HR) is the number of heartbeats per minute
- ◆ Heart Rate Variability (HRV) is the fluctuation in the time intervals between adjacent heartbeats
- ◆ HRV is concerned with analyzing the intervals between heart beats, which are called inter-beat intervals (IBIs)



Heart Rate Variability





Figure 2: Electrophysiology of the heart (redrawn from Malmivuo & Plonsey 1995). The different waveforms for each of the specialized cells found in the heart are shown. The latency shown approximates that normally found in the healthy heart.

- ◆ HRV is a complex pattern of oscillates (vibrations/frequencies)
- ◆ Examples: Heart Rate, Blood Pressure, Body Temperature, Energy Level, Mood, Circadian Rhythms
- ◆ Healthy biological systems exhibit complex patterns of variability that can be described by mathematical chaos
- ◆ The variability of non-linear systems provides the flexibility to rapidly cope with an uncertain and changing environment
- ◆ The "oscillation source" of the HRV signal is the time between R-waves or QRS complex
- ◆ Oscillations in one cardiovascular function causes identical oscillations in others

Heart Rate Variability



- ◆ HRV is conducting an ANS analysis through the portal of the heart
- ◆ HRV is a commonly used tool when trying to assess the function of cardiac autonomic regulation
- ◆ The ease at which an individual can transition between high and low arousal states is dependent on the ability of the ANS to rapidly vary heart rate
- ◆ More HRV implies more PNS activity (vagal tone) and better recovery or functional homeostatic capabilities
- ◆ And more homeostatic capacities implies better overall health and wellness

Exagon Sense

- ◆ Photoplethysmographic sensor (2-way: HR and SP02)
- ◆ Refined and very accurate algorithm (Histogram)
- ◆ Dynamic adjustment of applied magnetic field intensity (every 60-90 heart beats)
- ◆ Data Export, KUBIOS-compatible (www.kubios.com)
- ◆ Evaluation of HRV properties
- ◆ Exagon Sense in combination with KUBIOS can be utilized as a monitoring unit to observe and assess the autonomic nervous system over time



Heart Rate Variability

Heart rate variability



For scientific research and professional use


Kubios HRV Premium

- ✓ Supports IBI, ECG & PPG data
- ✓ Advanced pre-processing
- ✓ The most detailed HRV analysis on the market
- ✓ Detailed reports and exporting



Supports your measurement device

Kubios HRV Premium is compatible with several commonly used electrocardiogram (ECG) devices and with most of the HR monitors on the market. The only requirement for the HR monitor is that it can record inter-beat interval (IBI) data.




Data from your device

www.kubios.com

1 Niskanen J-P et al. Software for advanced HRV analysis. *Comp Meth Programs Biomed*, 76(1):73-81, 2004.
2 Tarvainen MP et al. Kubios HRV – Heart rate variability analysis software. *Comp Meth Programs Biomed*, 113(1):210-220, 2014.

- ◆ Kubios was founded in 1996, market leader in HRV analysis software for scientific research and professional use
- ◆ Kubios customers: Researchers, Well-being Therapists, Sports/Exercise Coaches, Athletes and Home Users for personal monitoring, etc.
- ◆ Applications of HRV include; Medical Research, Stress, Well-being, Fitness and Exercise
- ◆ Applications that relate to the iMRS Prime include; Stress Assessment, HRV Biofeedback, Resonant Frequencies, Physiotherapy and Yoga, etc.

Heart Rate Variability



- ◆ These graphs give a quick view about the level of HRV subjects with respect to normal values
- ◆ The interpretation of the PNS and SNS index is straightforward
- ◆ A PNS (or SNS) index value of zero means that the parameters reflecting parasympathetic (or sympathetic) activity are on average equal to normal population average